




# PRODUCT CATALOG

---

#1 DIGITAL IMPLANT COMPANY

# PEAK PERFORMANCE AT THE RIGHT PRICE.



## TRI® PERFORMANCE CONCEPT

The TRI® Performance Concept integrates 4 Gold Standards, in order to make the life of the practitioner as **SIMPLE** and yet **HIGH PERFORMING** and to achieve natural esthetic results in daily practice.

The TRI® Performance Concept is incorporated in the **Bone Level** (TRI®-Narrow, TRI®-Vent) and **Tissue Level** (TRI®-Octa) implant lines and ensures optimal esthetic results both in the anterior and posterior area.

Designed and Manufactured by TRI® in Switzerland 

# A LEAN AND INTELLIGENT IMPLANT CONCEPT.



The **TRI® Dental Implant System** combines the pure simplicity of one connection with the esthetic flexibility of the **TRI® Soft Tissue Concept**. A unique lean, and yet intelligent blend for success in daily practice.



## TRI®-Vent

TRI®-VENT IMPLANTS - Includes Surgical Cover Screw TVSCS								
Implant Ø	Implant Platform	Connection	6.5 mmL	8 mmL	10 mmL	11.5 mmL	13 mmL	16 mmL
3.75 mm	3.5 mmD	⬡ 2.5 mmD Internal Hexagon		TV37M08	TV37M10	TV37M11	TV37M13	TV37M16
4.1 mm			TV41M06	TV41M08	TV41M10	TV41M11	TV41M13	TV41M16
4.7 mm			TV47M06	TV47M08	TV47M10	TV47M11	TV47M13	TV47M16



## TRI®-Narrow

TRI®-NARROW IMPLANTS - Includes Surgical Cover Screw TNSCS							
Implant Ø	Implant Platform	Connection		10 mmL	11 mmL	13 mmL	16 mmL
3.3 mm	3.2 mmD	⬡ 2.1 mmD Internal Hexagon		TN33M10	TN33M11	TN33M13	TN33M16



## TRI®-Octa

TRI®-OCTA IMPLANTS - Includes Surgical Cover Screw TOSCS							
Implant Ø	Implant Platform	Connection	6.5 mmL	8 mmL	10 mmL	11.5 mmL	13 mmL
3.75 mm	4.8 mmD 1.8 mm Shoulder	⬡ 3.1 mmD Internal Octagon		TO37M08	TO37M10	TO37M11	TO37M13
4.1 mm			TO41M06	TO41M08	TO41M10	TO41M11	TO41M13
4.7 mm			TO47M06	TO47M08	TO47M10	TO47M11	TO47M13

# TRI®- VENT PROSTHETICS

## 3.5 mmD Implant Platform

### Healing Collars

#### TISSUE HEALING - HEALING COLLARS

Product Code		Ø	GH	Material	Qty	Sterile
TVHC35-30	Healing Collar	3.5 mm	3.0 mm	Ti-6Al-4V	1	Ⓡ
TVHC35-45	Healing Collar	3.5 mm	4.5 mm	Ti-6Al-4V	1	Ⓡ
TVHC35-60	Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1	Ⓡ
TVHC35-60-A	Suture Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1	Ⓡ
TVHC50-30	Healing Collar	5.0 mm	3.0 mm	Ti-6Al-4V	1	Ⓡ
TVHC50-45	Healing Collar	5.0 mm	4.5 mm	Ti-6Al-4V	1	Ⓡ
TVHC50-60	Healing Collar	5.0 mm	6.0 mm	Ti-6Al-4V	1	Ⓡ
TVHC50-60-A	Suture Healing Collar	5.0 mm	6.0 mm	Ti-6Al-4V	1	Ⓡ



#### SURGICAL COVER SCREWS

Product Code		Ø	GH	Material	Qty	Sterile
TVSCS	Replacement Cover Screw	3.5 mm		Ti-6Al-4V	1	Ⓡ
TVSCS-SL	Cover Screw for Sinus Lift	5.2 mm		Ti-6Al-4V	1	Ⓡ



### Impression Transfer Components, Implant Analogs

#### DIRECT TRANSFER FOR OPEN-TRAY IMPRESSIONS TRI® 3D-TOUCH IMPRESSION - INCLUDES SCREW DTRS-TV05

Product Code		Ø	GH	Material	Qty	Sterile
TV05-35	Direct Transfer Component	3.5 mm		Ti-6Al-4V	1	
TV05-50	Direct Transfer Component	5.0 mm		Ti-6Al-4V	1	
DTRS-TV05	Direct Transfer Retaining Screw			Ti-6Al-4V	1	



#### INDIRECT TRANSFER / STRAIGHT ABUTMENT FOR CLOSED-TRAY IMPRESSIONS TRI® 2IN1 IMPRESSION ABUTMENTS - INCLUDES SCREWS RS-TV10 AND RS-TV10-ITC

Product Code		Ø	GH	Material	Qty	Sterile
TV10-35-ITC	Indirect Transfer Component	3.5 mm		Ti-6Al-4V	1	
TV10-50-ITC	Indirect Transfer Component	5.0 mm		Ti-6Al-4V	1	
RS-TV10-ITC	Indirect Transfer Component Impression Screw			POM	5	
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	



#### IMPLANT ANALOG

Product Code		Material	Qty	Sterile
TVA-1	Digital Implant Analog	Steel	1	
TVA-5	Digital Implant Analog (5 Pack)	Steel	5	
TVIA-1	Implant Analog	Steel	1	
TVIA-5	Implant Analog (5 Pack)	Steel	5	



### Components for Cement-Retained Restorations


#### PLASTIC TEMPORARY ABUTMENTS - INCLUDES SCREW RS-TV20

Product Code		Ø	Material	Qty	Sterile
TV10-PTA	Temporary Abutment, Straight	4.5 mm	PEEK	1	
RS-TV20	Replacement Retaining Screw		Ti-6Al-4V	1	
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Red		Ti-6Al-4V	1	
TV10-EHA-50	Esthetic Healing Abutment	5.0 mm	PEEK	1	Ⓡ
TV10-EHA-70	Esthetic Healing Abutment	7.0 mm	PEEK	1	Ⓡ
RS-TV10	Replacement Retaining Screw		Ti-6Al-4V	1	
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green		Ti-6Al-4V	1	




#### STRAIGHT ABUTMENTS WITH TRI®-FRICTION - INCLUDES SCREW RS-TV10

Product Code		Ø	GH	Material	Qty	Sterile
TV10-35-F	Straight Abutment with TRI®-Friction	3.5 mm		Ti-6Al-4V	1	
TV10-50-F	Straight Abutment with TRI®-Friction	5.0 mm		Ti-6Al-4V	1	
TV10-41-F	Contour Abutment, Straight with TRI®-Friction	4.5 mm	1.0 mm	Ti-6Al-4V	1	
TV10-43-F	Contour Abutment, Straight with TRI®-Friction	4.5 mm	3.0 mm	Ti-6Al-4V	1	
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	



#### ANGLED ABUTMENTS WITH TRI®-FRICTION - INCLUDES SCREWS RS-TV10/RS-TV20

Product Code		Ø	GH	Material	Qty	Sterile
TV20-50-F	Angled Abutment, Standard - 20° with TRI®-Friction	5.0 mm		Ti-6Al-4V	1	
RS-TV20	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Red			Ti-6Al-4V	1	
TV20-41-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.5 mm	1.0 mm	Ti-6Al-4V	1	
TV20-43-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.5 mm	3.0 mm	Ti-6Al-4V	1	
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	




CH = Gingiva Height

Castable Abutments

CASTABLE ABUTMENTS, ENGAGING - INCLUDES SCREW RS-TV10

Product Code		Ø	GH	Material	Qty	Sterile
TV30-50-F	Gold Castable Abutment with TRI®-Friction	5.0 mm		HMA	1	
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	
TV35-50	Platin-Iridium Castable Abutment	5.0 mm		PtIr	1	
RS-TV20	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	



Components for Screw-Retained Restorations

SCREW-RETAINED MULTI UNIT ABUTMENTS

Product Code		Ø	GH	Material	Qty	Sterile
TV40-01	Screw-Retained Abutment, Straight	4.5 mm	1.0 mm	Ti-6Al-4V	1	
TV40-02	Screw-Retained Abutment, Straight	4.5 mm	2.0 mm	Ti-6Al-4V	1	
TV40-04	Screw-Retained Abutment, Straight	4.5 mm	4.0 mm	Ti-6Al-4V	1	
TV40-06	Screw-Retained Abutment, Straight	4.5 mm	6.0 mm	Ti-6Al-4V	1	
TV50-17	Screw Retained Abutment, 17° angled	4.5 mm	2.0 mm	Ti-6Al-4V	1	
TV50-30	Screw Retained Abutment 30° Angled	4.5 mm	3.0 mm	Ti-6Al-4V	1	
HC-SRAS	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1	
TV40-SCAN	Abutment Level Scanbody			PEEK	1	
TTA	Titanium Abutment incl. Screw RS-TTA			Ti-6Al-4V	1	
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1	
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1	
DTC-SRA	Direct Transfer Component for Screw-Retained Abutment incl. TRS-SRA			Ti-6Al-4V	1	
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1	
IA-SRA	Tapered Abutment Replica			Ti-6Al-4V	1	
MUA	Multi Unit Digital Implant Analog			Ti-6Al-4V	1	
MUA-5	Multi Unit Digital Implant Analog			Ti-6Al-4V	5	
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4	
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1	
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1	
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1	



Overdenture Restorations

BALL ABUTMENTS - INCLUDES CAHL

Product Code		Ø	GH	Material	Qty	Sterile
TV60-02	Ball Abutment incl. CAHL	4.5 mm	2.0 mm	Ti-6Al-4V	1	
TV60-04	Ball Abutment incl. CAHL	4.5 mm	4.0 mm	Ti-6Al-4V	1	
CAHL	Cap Attachment Housing and Liner - Transparent	4.0 mm		Stainless steel	1	
CAH	Cap Attachment Housing			Ti	1	
CALT	Cap Attachment Liner- (Transparent - Standard retention)			Nylon	1	
CALG	Cap Attachment Liner - (Gray - Rigid retention)			Nylon	1	
CALB	Cap Attachment Liner - Lab Use - Black			Nylon	1	



TRI® RAPTOR ABUTMENTS - INCLUDES CAPS ASSORTMENT KIT (PROD. CODE 192ECE)


Product Code		Ø	GH	Material	Qty	Sterile
130CV.01	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		1.0 mm	TiN	1	
130CV.02	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		2.0 mm	TiN	1	
130CV.03	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		3.0 mm	TiN	1	
130CV.04	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		4.0 mm	TiN	1	
130CV.05	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		5.0 mm	TiN	1	
130CV.06	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		6.0 mm	TiN	1	
192ECE	Caps Assortment Kit			Ti, Nylon	1	
335SBC	Smart Box Kit			Ti, Nylon	1	
144AE	Analog for Lab			Stainless	2	
044CAIN	Transfer for Impression			Stainless	2	



Components for TRI + Digital CAD/CAM

CAD/CAM ABUTMENTS - INCLUDES SCREW RS-TV10

Product Code		Ø	GH	Material	Qty	Sterile
TV70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1	
TV70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.5 mm	0.7 mm	Ti-6Al-4V	1	
TV70-20-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.5 mm	2.0 mm	Ti-6Al-4V	1	
TV70-07-NE	Non-Engaging Titatium Bonding-Base	4.5 mm	0.7 mm	Ti-6Al-4V	1	
TV70-20-NE	Non-Engaging Titatium Bonding-Base	4.5 mm	2.0 mm	Ti-6Al-4V	1	
TV70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1	
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	




# TRI®- NARROW PROSTHETICS

## 3.2 mmD Implant Platform

### Healing Collars

#### TISSUE HEALING - HEALING COLLARS

Product Code		Ø	GH	Material	Qty	Sterile
TNHC35-30	Healing Collar	3.5 mm	3.0 mm	Ti-6Al-4V	1	
TNHC35-45	Healing Collar	3.5 mm	4.5 mm	Ti-6Al-4V	1	
TNHC35-60	Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1	
TNHC35-60-A	Suture Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1	
TNHC45-30	Healing Collar	4.5 mm	3.0 mm	Ti-6Al-4V	1	
TNHC45-50	Healing Collar	4.5 mm	5.0 mm	Ti-6Al-4V	1	



#### SURGICAL COVER SCREW

Product Code		Ø	GH	Material	Menge	Sterile
TNSCS	Replacement Cover Screw	3.2 mm		Ti-6Al-4V	1	



### Impression Transfer Components, Implant Analogs

#### DIRECT TRANSFER FOR OPEN-TRAY IMPRESSIONS INCLUDES SCREW DTRS-TN05

Product Code		Ø	GH	Material	Qty	Sterile
TN05-35	Direct Transfer Component	3.5 mm		Ti-6Al-4V	1	
DTRS-TN05	Direct Transfer Retaining Screw			Ti-6Al-4V	1	



#### INDIRECT TRANSFER / STRAIGHT ABUTMENT FOR CLOSED-TRAY IMPRESSIONS TRI® 2IN1 IMPRESSION ABUTMENT – INCLUDES SCREW RS-TN10

Product Code		Ø	GH	Material	Qty	Sterile
TN10-35-ITC	Indirect Transfer Component	3.5 mm		Ti-6Al-4V	1	
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1	



#### IMPLANT ANALOG

Product Code		Material	Qty	Sterile
TNA-1	Digital Implant Analog	Ti-6Al-4V	1	
TNA-5	Digital Implant Analog (5 Pack)	Ti-6Al-4V	5	
TNIA-1	Implant Analog	Ti-6Al-4V	1	
TNIA-5	Implant Analog (5 Pack)	Ti-6Al-4V	5	



### Components for Cement-Retained Restorations


#### PLASTIC TEMPORARY ABUTMENT - INCLUDES SCREW RS-TN20

Product Code		Ø	GH	Material	Qty	Sterile
TN10-PTA	Temporary Abutment, Straight	3.5 mm		PEEK	1	
RS-TN20	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TN20-Lab	Replacement Retaining Screw - Lab Use - Red			Ti-6Al-4V	1	




#### STRAIGHT ABUTMENTS WITH TRI®-FRICTION - INCLUDES SCREW RS-TN10

Product Code		Ø	GH	Material	Qty	Sterile
TN10-35-F	Straight Abutment with TRI®-Friction	3.5 mm		Ti-6Al-4V	1	
TN10-41-F	Contour Abutment, Straight with TRI®-Friction	4.0 mm	1.0 mm	Ti-6Al-4V	1	
TN10-43-F	Contour Abutment, Straight with TRI®-Friction	4.0 mm	3.0 mm	Ti-6Al-4V	1	
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	



#### ANGLED ABUTMENTS WITH TRI®-FRICTION - INCLUDES SCREW RS-TN10/RS-TN20

Product Code		Ø	GH	Material	Qty	Sterile
TN20-45-F	Angled Abutment, 20° with TRI®-Friction	4.5 mm		Ti-6Al-4V	1	
TN20-41-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.0 mm	1.0 mm	Ti-6Al-4V	1	
TN20-43-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.0 mm	3.0 mm	Ti-6Al-4V	1	
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	



### Castable Abutments

#### GOLD ABUTMENTS, ENGAGING - INCLUDES SCREW RS-TN10

Product Code		Ø	GH	Material	Qty	Sterile
TN30-45-F	Gold Castable Abutment with TRI®-Friction	4.5 mm		HMA	1	
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	



### Components for TRI+ Digital CAD/CAM

#### CAD/CAM ABUTMENTS - INCLUDES SCREW RS-TN10




Product Code		Ø	GH	Material	Qty	Sterile
TN70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1	
TN70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.2 mm	0.7 mm	Ti-6Al-4V	1	
TN70-07-NE	Non-Engaging Titanium Bonding-Base	4.2 mm	0.7 mm	Ti-6Al-4V	1	
TN70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1	
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1	
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1	




GH = Gingiva Height

Components for Screw-Retained Restorations

SCREW-RETAINED MULTI UNIT ABUTMENTS

Product Code		Ø	GH	Material	Qty	Sterile
TN40-02	Screw-Retained Abutment, Straight	4.5 mm	2.0 mm	Ti-6Al-4V	1	
TN40-04	Screw-Retained Abutment, Straight	4.5 mm	4.0 mm	Ti-6Al-4V	1	
HC-SRAS	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1	
TV40-SCAN	Abutment Level Scanbody			PEEK	1	
TTA	Titanium Abutment incl. Screw RS-TTA			Ti-6Al-4V	1	
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1	
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1	
DTC-SRA	Direct Trasfer Component for Screw-Retained Abutment incl. DTRS			Ti-6Al-4V	1	
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1	
IA-SRA	Tapered Abutment Replica			Ti-6Al-4V	1	
MUA	Multi Unit Digital Implant Analog			Ti-6Al-4V	1	
MUA-5	Multi Unit Digital Implant Analog			Ti-6Al-4V	5	
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4	
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1	
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1	
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1	



Overdenture Restorations

BALL ABUTMENTS - INCLUDES CAHL

Product Code		Ø	GH	Material	Qty	Sterile
TN60-02	Ball Abutment incl. CAHL	4.5 mm	2.0 mm	Ti-6Al-4V	1	
TN60-04	Ball Abutment incl. CAHL	4.5 mm	4.0 mm	Ti-6Al-4V	1	
CAHL	Cap Attachment Housing and Liner - Transparent			Stainless steel	1	
CAH	Cap Attachment Housing			Stainless steel	1	
CALT	Cap Attachment Liner - Transparent			Nylon	1	
CALG	Cap Attachment Liner - Gray			Nylon	1	
CALB	Cap Attachment Liner - Lab Use - Black			Nylon	1	



TRI® RAPTOR ABUTMENTS - INCLUDES CAPS ASSORTMENT KIT (PROD. CODE 192ECE)




Product Code		Ø	GH	Material	Qty	Sterile
130TN.01	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		1.0 mm	TiN	1	
130TN.02	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		2.0 mm	TiN	1	
130TN.03	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		3.0 mm	TiN	1	
130TN.04	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		4.0 mm	TiN	1	
130TN.05	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		5.0 mm	TiN	1	
192ECE	Caps Assortment Kit			Ti, Nylon	1	
335SBC	Smart Box Kit			Ti, Nylon	1	
144AE	Analog for Lab			Stainless	2	
044CAIN	Transfer for Impression			Stainless	2	



TRI®- OCTA Prosthetics  
4.8 mmD Implant Platform

Healing Collars

TISSUE HEALING - HEALING COLLARS


Product Code		Ø	GH	Material	Qty	Sterile
TOSCS	Surgical Cover Screw	3.5 mm		Ti-6Al-4V	1	
TOHC51-15	Healing Collar	5.1 mm	1.5 mm	Ti-6Al-4V	1	
TOHC51-30	Healing Collar	5.1 mm	3.0 mm	Ti-6Al-4V	1	



Impression Transfer Components, Implant Analogs

DIRECT TRANSFER FOR OPEN-TRAY IMPRESSIONS  
TRI® 3D-TOUCH IMPRESSION - INCLUDES SCREW DTRS-TO05

Product Code		Ø	GH	Material	Qty	Sterile
TO05-51	Direct Transfer Component	5.1 mm		Ti-6Al-4V	1	
DTRS-TO05	Direct Transfer Retaining Screw			Ti-6Al-4V	1	



INDIRECT TRANSFER / STRAIGHT ABUTMENT FOR CLOSED-TRAY IMPRESSIONS  
TRI® 2IN1 IMPRESSION ABUTMENT - INCLUDES SCREW RS-TO

Product Code		Ø	GH	Material	Qty	Sterile
TO10-51-ITC	Indirect Transfer Component	5.1 mm		Ti-6Al-4V	1	
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1	



IMPLANT ANALOG

Product Code		Material	Qty	Sterile
TOA-1	Digital Implant Analog	Ti-6Al-4V	1	
TOA-5	Digital Implant Analog (5 Pack)	Ti-6Al-AV	5	
TOIA	Implant Analog	Ti-6Al-AV	1	
TOIA-5	Implant Analog (5 Pack)	Ti-6Al-AV	5	



GH = Gingiva Height

Components for Cement-Retained Restorations

PLASTIC TEMPORARY ABUTMENTS - INCLUDES SCREW RS-TO

Product Code		Ø	GH	Material	Qty	Sterile
TO10-PTA	Temporary Abutment, Straight	5.1 mm		PEEK	1	
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1	



TO10-PTA

STRAIGHT ABUTMENTS WITH TRI® FRICTION - INCLUDES SCREW RS-TO

Product Code		Ø	GH	Material	Qty	Sterile
TO10-35-F	Straight Abutment with TRI® Friction	3.5 mm		Ti-6Al-4V	1	
TO10-51-F	Straight Abutment with TRI® Friction	5.1 mm		Ti-6Al-4V	1	
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1	



TO10-35-F



TO10-51-F

ANGLED ABUTMENTS WITH TRI® FRICTION - INCLUDES SCREW RS-TO

Product Code		Ø	GH	Material	Qty	Sterile
TO20-51-F	Angled Abutment with TRI® Friction, Standard - 20°	5.1 mm		Ti-6Al-4V	1	
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1	



TO20-51-F

Components for Screw-Retained Restorations

SCREW-RETAINED MULTI UNIT ABUTMENTS

Product Code		Ø	GH	Material	Qty	Sterile
TO40-16	Screw-Retained Abutment, Straight	4.5 mm	1.6 mm	Ti-6Al-4V	1	Ⓡ
TO40-30	Screw-Retained Abutment, Straight	4.5 mm	3.0 mm	Ti-6Al-4V	1	Ⓡ
HC-SRAS	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V		Ⓡ
TV40-Scan	Abutment Level Scanbody			PEEK	1	
TTA	Titanium Abutment incl. Screw RS-TTA			Ti-6Al-4V	1	
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1	
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1	
DTC-SRA	Direct Trasfer Component for Screw-Retained Abutment incl. DTRS			Ti-6Al-4V	1	
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1	
IA-SRA	Tapered Abutment Replica			Ti-6Al-4V	1	
MUA	Multi Unit Digital Implant Analog			Ti-6Al-4V	1	
MUA-5	Multi Unit Digital Implant Analog			Ti-6Al-4V	5	
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4	
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1	
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1	
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1	



TO40-16



TO40-30



TTA



PCC



DTC-SRA



DTRS-SRA



IA-SRA

Overdenture Restorations

BALL ABUTMENTS - INCLUDES CAHL

Product Code		Ø	GH	Material	Qty	Sterile
TO60-16	Ball Abutment incl. CAHL	4.5 mm	1.6 mm	Ti-6Al-4V	1	
TO60-30	Ball Abutment incl. CAHL	4.5 mm	3.0 mm	Ti-6Al-4V	1	



TO60-16



TO60-30

CAHL	Cap Attachment Housing and Liner - Transparent			Stainless steel	1	
CAH	Cap Attachment Housing			Stainless steel	1	
CALT	Cap Attachment Liner - Transparent			Nylon	1	
CALG	Cap Attachment Liner - Gray			Nylon	1	
CALB	Cap Attachment Liner - Lab Use - Black			Nylon	1	



CAHL

TRI® RAPTOR - INCLUDES CAPS ASSORTMENT KIT (PROD. CODE 192ECE)

Product Code		Ø	GH	Material	Qty	Sterile
130BFT.00	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		0.0 mm	TiN	1	
130BFT.01	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		1.0 mm	TiN	1	
130BFT.02	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		2.0 mm	TiN	1	
130BFT.03	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		3.0 mm	TiN	1	
130BFT.04	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		4.0 mm	TiN	1	
130BFT.05	TRI® Raptor Abutment - includes Caps Assortment Kit (192ECE)		5.0 mm	TiN	1	
192ECE	Caps Assortment Kit			Ti, Nylon	1	
335SBC	Smart Box Kit			Ti, Nylon	1	
144AE	Analog for Lab			Stainless	2	
044CAIN	Transfer for Impression			Stainless	2	



192ECE



144AE



044CAIN

Components for TRI+ Digital CAD/CAM

CAD/CAM ABUTMENT - INCLUDES SCREW RS-TO

Product Code		Ø	GH	Material	Qty	Sterile
TO70-Scan	Scanbody for Laboratory and Intraoral Scanner			PEEK	1	
TO70-07-F	Engaging Titanium Bonding-Base with TRI® Friction	5.0 mm	0.7 mm	Ti-6Al-4V	1	
TO70-07-NE	Non-Engaging Titanium Bonding-Base	5.0 mm		Ti-6Al-4V	1	
TO70-Blank-F	Titanium Milling Blank with TRI®-Friction	5.0 mm		Ti-6Al-4V	1	
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1	



TO70-SCAN



TO70-07-F



TO70-07-NE





























TO70-Blank-F









GH = Gingiva Height

# INSTRUMENTS

## ONE SURGICAL KIT FOR ALL IMPLANT LINES.

### Surgical Kit for the TRI® Dental Implant System

Starting Drills		Description	
LD	Nr. 1	Lance Drill	
RB-2.3	Nr. 2	Round Bur	
DE	Nr. 3	Drill Extender	
Surgical Drills		Description	
TPD-2.3S	Nr. 4	Pilot Drill - Short - 2.3 mm	
TPD-2.3L	Nr. 5	Pilot Drill - Long - 2.3 mm	
TND-3.3L	Nr. 6	TRI®-Narrow Final Drill - Long - 2.4/2.8 mm	
TAP3.3	Nr. 7	TRI®-Narrow TAP	
SBD-3.7S	Nr. 8	3.7 Soft Bone Drill - Short - 2.7/3.1 mmD	
SBD-3.7L	Nr. 9	3.7 Soft Bone Drill - Long - 2.7/3.1 mmD	
DBD-3.7L	Nr. 10	3.7 Dense Bone Drill - Long - 2.9/3.6 mmD	
TAP3.7	Nr. 11	3.7 Tap	
SBD-4.1S	Nr. 12	4.1 Soft Bone Drill - Short - 3.0/3.4 mmD	
SBD-4.1L	Nr. 13	4.1 Soft Bone Drill - Long - 3.0/3.4 mmD	
DBD-4.1L	Nr. 14	4.1 Dense Bone Drill - Long – 3.2/3.9 mmD	
TAP4.1	Nr. 15	4.1 Tap	
SBD-4.7S	Nr. 16	4.7 Soft Bone Drill - Short - 3.6/4.0 mmD	
SBD-4.7L	Nr. 17	4.7 Soft Bone Drill - Long - 3.6/4.0 mmD	
DBD-4.7L	Nr. 18	4.7 Dense Bone Drill - Long – 3.8/4.6 mmD	
TAP4.7	Nr. 19	4.7 Tap	
PP-L	Nr. 21	Paralell Pin - Long	
Implant Drivers		Description	
VD -S		TRI®-Vent Implant Driver - Short	
VD-L		TRI®-Vent Implant Driver - Long	
ND-S		TRI®-Narrow Implant Driver - Short	
ND-L		TRI®-Narrow Implant Driver - Long	
OD-S		TRI®-Octa Implant Driver - Short	
OD-L		TRI®-Octa Implant Driver - Long	

Prosthetic Tools		
PD-ISO-S*	TRI® Prosthetic Driver for Motor - Short	
PD-ISO-L-Ang*	TRI® Prosthetic Driver Angulated - Long	
487,ICE	TRI® Raptor Insertion and Extractor Tool for caps	
ITFC	Insertion Tool for Caps-Ball Abutment	
RTFC	Removal Tool for Caps-Ball Abutment	
Torque Wrench		
TW	Torque Wrench	
TW - Adapter	Torque Wrench Adapter	
Drill Stop*		
TDS6	Drill Stop / 6 mm	
TDS8	Drill Stop / 8 mm	
TDS10	Drill Stop / 10 mm	
TDS11	Drill Stop / 11 mm	
TDS13	Drill Stop / 13 mm	
TPDS	Guided Surgery Sleeve for TPD-2.3L	
Other		
TSK	TRI® Surgical Kit - Complete	
TSK-Medium	TRI® Surgical Kit - Medium	
TSK-Tray	TRI® Surgical Kit - Tray	
TRI Sinus Kit	TRI® Sinus Lift Kit - Complete	
TRI Sinus Suction	TRI® Sinus Suction Device	



Prosthetic Kit for the TRI®Dental Implant System		
Catalog Number		
PK	TRI® Prosthetic Kit (includes: PD-ISO-L, PD-L, PD-S, TORT, TRT, TW, TW-Adapter)	
PK-Tray	Prosthetic Kit-Tray	
PD-ISO-L*	TRI®Prosthetic Driver for Motor-Long	
PD-L	TRI®Prosthetic Driver-Long	
PD-S	TRI®Prosthetic Driver-Short	
TORT	Removal Tool for TRI®-Octa Friction Abutments	
TRT	Removal Tool for TRI®-Vent & Narrow Friction Abutments	
TW	Torque Wrench	
TW-Adapter	Torque Wrench Adapter	



\* not included in TSK

# DRILL SEQUENCE

### 3.3 mmD TRI®-Narrow

Lance Drill

TPD-2.3L  
PilotDrill 2.3 mmD  
optional with  
Drill Stop

SOFT / DENSE BONE  
TND-3.3L  
2.4/2.8 mmD Drill

DENSE BONE  
TAP-3.3  
Cortical Bone TAP  
in D1 Bone

### 3.75 mmD TRI®-Vent and TRI®-Octa

Lance Drill

TPD-2.3L  
PilotDrill 2.3 mmD  
optional with  
Drill Stop

SOFT BONE  
SBD-3.7L  
2.7/3.1 mmD Drill

DENSE BONE  
DBD-3.7L  
2.9/3.6 mmD Drill

DENSE BONE  
TAP-3.7  
Cortical Bone TAP  
in D1 Bone

### 4.1 mmD TRI®-Vent and TRI®-Octa

Lance Drill

TPD-2.3L  
PilotDrill 2.3 mmD  
optional with  
Drill Stop

SBD-3.7L  
2.7/3.1 mmD Drill

SOFT BONE  
SBD-4.1L  
3.0/3.4 mmD Drill

DENSE BONE  
DBD-4.1L  
3.2/3.9 mmD Drill

DENSE BONE  
TAP-4.1  
Cortical Bone TAP  
in D1 Bone

### 4.7 mmD TRI®-Vent and TRI®-Octa

Lance Drill

TPD-2.3L  
PilotDrill 2.3 mmD  
optional with  
Drill Stop

SBD-3.7L  
2.7/3.1 mmD Drill


SBD-4.1L  
3.0/3.4 mmD Drill

SOFT BONE  
SBD-4.7L  
3.6/4.0 mmD Drill


DENSE BONE  
DBD-4.7L  
3.8/4.5 mmD Drill

DENSE BONE  
TAP-4.7  
Cortical Bone TAP  
in D1 Bone

# DRILL AND IMPLANT LENGTHS



# TORQUE RECOMMENDED



# TRI® POD




## PACKAGING HANDLING



**Step 1**  
Open Lid



**Step 2**  
Rotate To Remove Implant

- 3 Options for Implant Insertion**
-  **Ratchet**
  -  **Handpiece**
  -  **By hand**



**Step 3**  
Remove Cover Screw

# IMPLANTOLOGY

## MOTOR SYSTEM MD11

Unsurpassed reliability - the surgery motor system from 

The new MD 11 together with its brother, the MD 30 Motor System, are developed with the intention to put extra smoothness into your hands. The further development of its electronic motor and a sophisticated motor control are key players in this new generation of Nougavag motor systems. This reflects also the new design with smooth edges for easy cleaning and a display that gives you all the information at a glance. Single function keys with haptic and audible feedback guarantee for accurate device setup.



To learn more about our products visit our TRI® Academy Courses and Webinars: [academy@tri-implants.swiss](mailto:academy@tri-implants.swiss)

### Designed and Manufactured by TRI® in Switzerland

Designed in the heart of Switzerland, all implants are produced in specialized manufacturing units for dental implants. Therefore benefiting from more than 25 years of medical technology production know-how - the only way to guarantee Swiss Precision with high quality. Visit us in Switzerland: [visit@tri-implants.swiss](mailto:visit@tri-implants.swiss)





**TRI Dental Implants Int. AG**  
Switzerland  
Infoline: 00800 3313 3313  
[www.tri.swiss](http://www.tri.swiss)

Become a TRI® Follower



**TRI® Dental Implants**, the **TRI® Dental Implants** logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of the **TRI® Dental Implants Int. AG**. | 2020-02 | TRI® Product Catalog EN 2020-02/02



CE 0297 | FDA Registered