





**PRO-
DUCT**
CLASSIC
LINE **CATA-
LOG**



Narrow, Bone- and Tissue-Level

The Classic Line is incorporated in the Bone-Level (TRI[®]-Narrow, TRI[®]-Vent) and Tissue-Level (TRI[®]-Octa) implant lines and ensures optimal esthetic results both in the anterior and posterior area.

Lean & intelligent implant concept.



The TRI[®] Classic Line combines the pure simplicity of one connection with the esthetic flexibility of the TRI[®] Soft Tissue Concept. A unique lean, and yet intelligent blend for success in daily practice.



Surgical Kit

The TRI[®] Classic Line is characterized by one unique surgical kit for bone- and tissue-level line. Its intuitive and elegant design prioritizes efficiency and ease of use.

TRI[®]-Vent Implants - Includes Surgical Cover Screw TVSCS

Implant Ø	Platform	Connection	6.5mmL	8mmL	10mmL	11.5mmL	13mmL	16mmL
3.75 mm	3.5 mmD	2.5 mmD Internal Hexagon		TV37Mo8	TV37M10	TV37M11	TV37M13	TV37M16
4.1 mm			TV41Mo6	TV41Mo8	TV41M10	TV41M11	TV41M13	TV41M16
4.7 mm			TV47Mo6	TV47Mo8	TV47M10	TV47M11	TV47M13	TV47M16

TRI[®]-Narrow Implants - Includes Surgical Cover Screw TNSCS

Implant Ø	Platform	Connection	10mmL	11mmL	13mmL	16mmL
3.3 mm	3.0 mmD	2.1 mmD Internal Hexagon	TN33M10	TN33M11	TN33M13	TN33M16

TRI[®]-Octa Implants - Includes Surgical Cover Screw TOSCS



Implant Ø	Platform	Connection	6.5mmL	8mmL	10mmL	11.5mmL	13mmL
3.75 mm	4.8 mmD 1.8mm Shoulder	3.1 mmD Internal Octagon		TO37Mo8	TO37M10	TO37M11	TO37M13
4.1 mm			TO41Mo6	TO41Mo8	TO41M10	TO41M11	TO41M13
4.7 mm			TO47Mo6	TO47Mo8	TO47M10	TO47M11	TO47M13

TRI[®]-Vent Prosthetics


3.5 mmD Implant Platform

HEALING COLLARS

Tissue healing - Healing Collars, sterile



	Product Code		∅	GH	Material	Qty
	TVHC35-30	Healing Collar	3.5 mm	3.0 mm	Ti-6Al-4V	1
	TVHC35-45	Healing Collar	3.5 mm	4.5 mm	Ti-6Al-4V	1
	TVHC35-60	Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1
	TVHC50-30	Healing Collar	5.0 mm	3.0 mm	Ti-6Al-4V	1
	TVHC50-45	Healing Collar	5.0 mm	4.5 mm	Ti-6Al-4V	1
	TVHC50-60	Healing Collar	5.0 mm	6.0 mm	Ti-6Al-4V	1

Surgical cover screws - sterile

	Product Code		∅	GH	Material	Qty
	TVSCS	Replacement Cover Screw	3.5 mm		Ti-6Al-4V	1
	TVSCS-SL	Cover Screw for Sinus Lift	5.2 mm		Ti-6Al-4V	1



IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions TRI[®] 3D-Touch Impression - includes screw DTRS-TV05

	Product Code		∅	GH	Material	Qty
	TV05-35	Direct Transfer Component	3.5 mm		Ti-6Al-4V	1
	TV05-50	Direct Transfer Component	5.0 mm		Ti-6Al-4V	1
	DTRS-TV05	Direct Transfer Retaining Screw			Ti-6Al-4V	1

Indirect transfer/straight abutment for closed-tray impressions

TRI[®] 2in1 Impression Abutments - includes screws RS-TV10 and RS-TV10-ITC

	Product Code		∅	GH	Material	Qty
	TV10-35-ITC	Indirect Transfer Component	3.5 mm		Ti-6Al-4V	1
	TV10-50-ITC	Indirect Transfer Component	5.0 mm		Ti-6Al-4V	1
	RS-TV10-ITC	Indirect Transfer Component Impression Screw			POM	5
	RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
	RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

GH = Gingiva Height

Implant Analog



Product Code		Material	Qty
TVA-1	Digital Implant Analog	Steel	1
TVA-5	Digital Implant Analog (5 Pack)	Steel	5

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TV10/RS-TV20



Product Code		∅	Material	Qty
TV10-PTA	Temporary Abutment, Straight	4.5 mm	PEEK	1
RS-TV20	Replacement Retaining Screw		Ti-6Al-4V	1
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Red		Ti-6Al-4V	1
TV10-EHA-50*	Esthetic Healing Abutment	5.0 mm	PEEK	1
TV10-EHA-70*	Esthetic Healing Abutment	7.0 mm	PEEK	1
RS-TV10	Replacement Retaining Screw		Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green		Ti-6Al-4V	1

Straight Abutments with TRI[®]-Friction - includes screw RS-TV10



Product Code		∅	GH	Material	Qty
TV10-35-F	Straight Abutment with TRI [®] -Friction	3.5 mm		Ti-6Al-4V	1
TV10-50-F	Straight Abutment with TRI [®] -Friction	5.0 mm		Ti-6Al-4V	1
TV10-41-F	Contour Abutment, Straight with TRI [®] -Friction	4.5 mm	1.0 mm	Ti-6Al-4V	1
TV10-43-F	Contour Abutment, Straight with TRI [®] -Friction	4.5 mm	3.0 mm	Ti-6Al-4V	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

Angled Abutments with TRI[®]-Friction - includes screws RS-TV10/RS-TV20




Product Code		∅	GH	Material	Qty
TV20-50-F	Angled Abutment, Standard - 20° with TRI [®] -Friction	5.0 mm		Ti-6Al-4V	1
RS-TV20	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Red			Ti-6Al-4V	1
TV20-41-F	Contour Abutment, Angled - 17° with TRI [®] -Friction	4.5 mm	1.0 mm	Ti-6Al-4V	1
TV20-43-F	Contour Abutment, Angled - 17° with TRI [®] -Friction	4.5 mm	3.0 mm	Ti-6Al-4V	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

*sterile

CASTABLE ABUTMENTS

Castable Abutments, Engaging - includes screw RS-TV10/RS-TV20



Product Code		Ø	GH	Material	Qty
TV30-50-F	Gold Castable Abutment with TRI®-Friction	5.0 mm		HMA	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1
TV35-50	Platin-Iridium Castable Abutment	5.0 mm		PtIr	1
RS-TV20	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV20-Lab	Replacement Retaining Screw - Lab Use - Red			Ti-6Al-4V	1

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments




Product Code		Ø	GH	Material	Qty
TV40-01*	Screw-Retained Abutment, Straight	4.5 mm	1.0 mm	Ti-6Al-4V	1
TV40-02*	Screw-Retained Abutment, Straight	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV40-04*	Screw-Retained Abutment, Straight	4.5 mm	4.0 mm	Ti-6Al-4V	1
TV40-06*	Screw-Retained Abutment, Straight	4.5 mm	6.0 mm	Ti-6Al-4V	1
TV50-17*	Screw Retained Abutment, 17° Angled	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV50-30*	Screw Retained Abutment, 30° Angled	4.5 mm	3.0 mm	Ti-6Al-4V	1
HC-SRAS*	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1
TV40-SCAN	Abutment Level Scanbody			PEEK	1
TTA	Titanium Temporary Abutment incl. Screw RS-TTA			Ti-6Al-4V	1
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1
DTC-SRA	Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA			Ti-6Al-4V	1
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1
MUA	Multi-Unit Digital Implant Analog			Ti-6Al-4V	1
MUA-5	Multi-Unit Digital Implant Analog			Ti-6Al-4V	5
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1

*sterile

OVERDENTURE RESTORATIONS

Ball Abutments - includes CAHL



Product Code		Ø	GH	Material	Qty
TV60-02	Ball Abutment incl. CAHL	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV60-04	Ball Abutment incl. CAHL	4.5 mm	4.0 mm	Ti-6Al-4V	1
CAHL	Cap Attachment Housing and Liner - Transparent			Stainless steel	1
CAH	Cap Attachment Housing			Stainless steel	1
CALT	Cap Attachment Liner - (Transparent - Standard retention)			Nylon	1
CALG	Cap Attachment Liner - (Gray - Rigid retention)			Nylon	1
CALB	Cap Attachment Liner - Lab Use - Black			Nylon	1

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)



Product Code		Ø	GH	Material	Qty
130CV.01	TRI® Raptor Abutment		1.0 mm	TiN	1
130CV.02	TRI® Raptor Abutment		2.0 mm	TiN	1
130CV.03	TRI® Raptor Abutment		3.0 mm	TiN	1
130CV.04	TRI® Raptor Abutment		4.0 mm	TiN	1
130CV.05	TRI® Raptor Abutment		5.0 mm	TiN	1
130CV.06	TRI® Raptor Abutment		6.0 mm	TiN	1
192ECE	Caps Assortment Kit			Ti, Nylon	1
335SBC	Smart Box Kit			Ti, Nylon	1
144AE	Analog for Lab			Stainless	2
044CAIN	Transfer for Impression			Stainless	2

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TV10



Product Code		Ø	GH	Material	Qty
TV70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1
TV70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.5 mm	0.7 mm	Ti-6Al-4V	1
TV70-20-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV70-07-NE	Non-Engaging Titanium Bonding-Base	4.5 mm	0.7 mm	Ti-6Al-4V	1
TV70-20-NE	Non-Engaging Titanium Bonding-Base	4.5 mm	2.0 mm	Ti-6Al-4V	1
TV70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1
RS-TV10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TV10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

GH = Gingiva Height

TRI®-Narrow Prosthetics

3.0 mmD Implant Platform

HEALING COLLARS

Tissue healing - Healing Collars, sterile



TNHC35-30

Product Code		Ø	GH	Material	Qty
TNHC35-30	Healing Collar	3.5 mm	3.0 mm	Ti-6Al-4V	1
TNHC35-45	Healing Collar	3.5 mm	4.5 mm	Ti-6Al-4V	1
TNHC35-60	Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1
TNHC45-30	Healing Collar	4.5 mm	3.0 mm	Ti-6Al-4V	1
TNHC45-50	Healing Collar	4.5 mm	5.0 mm	Ti-6Al-4V	1

Surgical cover screws, sterile



TNCS

Product Code		Ø	GH	Material	Qty
TNCS	Replacement Cover Screw	3.2 mm		Ti-6Al-4V	1

IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions - TRI® 3D-Touch Impression - includes screw DTRS-TN05



TN05-35



DTRS-TN05

Product Code		Ø	GH	Material	Qty
TN05-35	Direct Transfer Component	3.5 mm		Ti-6Al-4V	1
DTRS-TN05	Direct Transfer Retaining Screw			Ti-6Al-4V	1

Indirect transfer/straight abutment for closed-tray impressions

TRI® 2in1 Impression Abutments - includes screw RS-TN10



TN 10-35-ITC

Product Code		Ø	GH	Material	Qty
TN10-35-ITC	Indirect Transfer Component	3.5 mm		Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1

Implant Analog



TNA



TNIA

Product Code		Material	Qty
TNA-1	Digital Implant Analog	Ti-6Al-4V	1
TNA-5	Digital Implant Analog (5 Pack)	Ti-6Al-4V	5

GH = Gingiva Height

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TN20



TN10-PTA

Product Code		Ø	Material	Qty
TN10-PTA	Temporary Abutment, Straight	3.5 mm	PEEK	1
RS-TN20	Replacement Retaining Screw		Ti-6Al-4V	1
RS-TN20-Lab	Replacement Retaining Screw - Lab Use - Red		Ti-6Al-4V	1

Straight Abutments with TRI®-Friction - includes screw RS-TN10



TN10-41-F



TN10-35-F

Product Code		Ø	GH	Material	Qty
TN10-35-F	Straight Abutment with TRI®-Friction	3.5 mm		Ti-6Al-4V	1
TN10-41-F	Contour Abutment, Straight with TRI®-Friction	4.0 mm	1.0 mm	Ti-6Al-4V	1
TN10-43-F	Contour Abutment, Straight with TRI®-Friction	4.0 mm	3.0 mm	Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

Angled Abutments with TRI®-Friction - includes screws RS-TN10



TN20-41-F



TN20-45-F

Product Code		Ø	GH	Material	Qty
TN20-45-F	Angled Abutment, Standard - 20° with TRI®-Friction	4.5 mm		Ti-6Al-4V	1
TN20-41-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.0 mm	1.0 mm	Ti-6Al-4V	1
TN20-43-F	Contour Abutment, Angled - 17° with TRI®-Friction	4.0 mm	3.0 mm	Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

CASTABLE ABUTMENTS

Gold abutments, engaging - includes screw RS-TN10



TN30-45-F

Product Code		Ø	Material	Qty
TN30-45-F	Gold Castable Abutment with TRI®-Friction	4.5 mm	HMA	1
RS-TN10	Replacement Retaining Screw		Ti-6Al-4V	1
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green		Ti-6Al-4V	1

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TN10



TN70-SCAN



TN70-07-F



TN70-07-NE



TN70-Blank-F

Product Code		Ø	GH	Material	Qty
TN70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1
TN70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	4.2 mm	0.7 mm	Ti-6Al-4V	1
TN70-07-NE	Non-Engaging Titanium Bonding-Base	4.2 mm	0.7 mm	Ti-6Al-4V	1
TN70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1
RS-TN10	Replacement Retaining Screw			Ti-6Al-4V	1
RS-TN10-Lab	Replacement Retaining Screw - Lab Use - Green			Ti-6Al-4V	1

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments



Product Code		Ø	GH	Material	Qty
TN40-02*	Screw-Retained Abutment, Straight	4.5 mm	2.0 mm	Ti-6Al-4V	1
TN40-04*	Screw-Retained Abutment, Straight	4.5 mm	4.0 mm	Ti-6Al-4V	1
HC-SRAS*	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1
TV40-SCAN	Abutment Level Scanbody			PEEK	1
TTA	Titanium Abutment incl. Screw RS-TTA			Ti-6Al-4V	1
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1
DTC-SRA	Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA			Ti-6Al-4V	1
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1
MUA	Multi-Unit Digital Implant Analog			Ti-6Al-4V	1
MUA-5	Multi-Unit Digital Implant Analog			Ti-6Al-4V	5
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1

OVERDENTURE RESTORATIONS

Ball Abutments - includes CAHL



Product Code		Ø	GH	Material	Qty
TN60-02	Ball Abutment incl. CAHL	4.5 mm	2.0 mm	Ti-6Al-4V	1
TN60-04	Ball Abutment incl. CAHL	4.5 mm	4.0 mm	Ti-6Al-4V	1
CAHL	Cap Attachment Housing and Liner - Transparent			Stainless steel	1
CAH	Cap Attachment Housing			Stainless steel	1
CALT	Cap Attachment Liner - (Transparent - Standard retention)			Nylon	1
CALG	Cap Attachment Liner - (Gray - Rigid retention)			Nylon	1
CALB	Cap Attachment Liner - Lab Use - Black			Nylon	1

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)



Product Code		Ø	GH	Material	Qty
130TN.01	TRI® Raptor Abutment		1.0 mm	TiN	1
130TN.02	TRI® Raptor Abutment		2.0 mm	TiN	1
130TN.03	TRI® Raptor Abutment		3.0 mm	TiN	1
130TN.04	TRI® Raptor Abutment		4.0 mm	TiN	1
130TN.05	TRI® Raptor Abutment		5.0 mm	TiN	1
192ECE	Caps Assortment Kit			Ti, Nylon	1
335SBC	Smart Box Kit			Ti, Nylon	1
144AE	Analog for Lab			Stainless	2
044CAIN	Transfer for Impression			Stainless	2

*sterile

TRI®-Octa Prosthetics 4.8 mmD Implant Platform

HEALING COLLARS


Tissue healing - Healing Collars, sterile



Product Code		Ø	GH	Material	Qty
TOSCS	Surgical Cover Screw	3.5 mm		Ti-6Al-4V	1
TOHC51-15	Healing Collar	5.1 mm	1.5 mm	Ti-6Al-4V	1
TOHC51-30	Healing Collar	5.1 mm	3.0 mm	Ti-6Al-4V	1

IMPRESSION TRANSFER COMPONENTS, IMPLANT ANALOGS

Direct transfer for open-tray impressions - TRI® 3D-Touch Impression - includes screw DTRS-TO05



Product Code		Ø	GH	Material	Qty
TO05-51	Direct Transfer Component	5.1 mm		Ti-6Al-4V	1
DTRS-TO05	Direct Transfer Retaining Screw			Ti-6Al-4V	1

Indirect transfer/straight abutment for closed-tray impressions

TRI® 2in1 Impression Abutments - includes screw RS-TO



TO10-51-ITC

Product Code		Ø	GH	Material	Qty
TO10-51-ITC	Indirect Transfer Component	5.1 mm		Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

Implant Analog



Product Code		Material	Qty
TOA-1	Digital Implant Analog	Ti-6Al-4V	1
TOA-5	Digital Implant Analog (5 Pack)	Ti-6Al-4V	5

GH = Gingiva Height

COMPONENTS FOR CEMENT-RETAINED RESTORATIONS

Plastic Temporary Abutments - includes screw RS-TO



TO10-PTA

Product Code		Ø	Material	Qty
TO10-PTA	Temporary Abutment, Straight	5.1 mm	PEEK	1
RS-TO	Replacement Retaining Screw		Ti-6Al-4V	1

Straight Abutments with TRI® Friction - includes screw RS-TO



TO10-35-F



TO10-51-F

Product Code		Ø	GH	Material	Qty
TO10-35-F	Straight Abutment with TRI®-Friction	3.5 mm		Ti-6Al-4V	1
TO10-51-F	Straight Abutment with TRI®-Friction	5.1 mm		Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

Angled Abutments with TRI®-Friction - includes screws RS-TO



TO20-51-F

Product Code		Ø	GH	Material	Qty
TO20-51-F	Angled Abutment, Standard - 20° with TRI®-Friction	5.1 mm		Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

COMPONENTS FOR SCREW-RETAINED RESTORATIONS

Screw-retained Multi-Unit Abutments



TO40-16



TO40-30



TTA



MUA



PCC



DTC-SRA



DTRS-SRA

Product Code		Ø	GH	Material	Qty
TO40-16*	Screw-Retained Abutment, Straight	4.5 mm	1.6 mm	Ti-6Al-4V	1
TO40-30*	Screw-Retained Abutment, Straight	4.5 mm	3.0 mm	Ti-6Al-4V	1
HC-SRAS*	Healing Cap for Screw-Retained Abutment			Ti-6Al-4V	1
TV40-SCAN	Abutment Level Scanbody			PEEK	1
TTA	Titanium Abutment incl. Screw RS-TTA			Ti-6Al-4V	1
RS-TTA	Replacement Retaining Screw for TTA			Ti-6Al-4V	1
RS-TTAL	Replacement Retaining Screw Long for TTA			Ti-6Al-4V	1
DTC-SRA	Direct Transfer Component for Screw-Retained Abutment incl. DTRS-SRA			Ti-6Al-4V	1
DTRS-SRA	Direct Transfer Retaining Screw			Ti-6Al-4V	1
MUA	Multi-Unit Digital Implant Analog			Ti-6Al-4V	1
MUA-5	Multi-Unit Digital Implant Analog			Ti-6Al-4V	5
RS-MUA-D-4	Replacement Retaining Screw for MUA			Ti-6Al-4V	4
PCC	Plastic Castable Coping incl. Screw RS-PCC			POM	1
RS-PCC	Replacement Retaining Screw for PCC			Ti-6Al-4V	1
RS-PCC-Lab	Replacement Retaining Screw - Lab Use - Blue			Ti-6Al-4V	1

OVERDENTURE RESTORATIONS

Ball Abutments - includes CAHL



TO60-16



TO60-30



CAHL

Product Code		Ø	GH	Material	Qty
To60-16	Ball Abutment incl. CAHL	4.5 mm	1.6 mm	Ti-6Al-4V	1
TO60-30	Ball Abutment incl. CAHL	4.5 mm	3.0 mm	Ti-6Al-4V	1
CAHL	Cap Attachment Housing and Liner - Transparent			Stainless steel	1
CAH	Cap Attachment Housing			Stainless steel	1
CALT	Cap Attachment Liner - (Transparent - Standard retention)			Nylon	1
CALG	Cap Attachment Liner - (Gray - Rigid retention)			Nylon	1
CALB	Cap Attachment Liner - Lab Use - Black			Nylon	1

TRI® Raptor Abutments - includes Caps Assortment Kit (Prod. Code 192ECE)



192ECE



144AE



044CAIN

Product Code		Ø	GH	Material	Qty
130BFT.00	TRI® Raptor Abutment		0.0 mm	TiN	1
130BFT.01	TRI® Raptor Abutment		1.0 mm	TiN	1
130BFT.02	TRI® Raptor Abutment		2.0 mm	TiN	1
130BFT.03	TRI® Raptor Abutment		3.0 mm	TiN	1
130BFT.04	TRI® Raptor Abutment		4.0 mm	TiN	1
130BFT.05	TRI® Raptor Abutment		5.0 mm	TiN	1
192ECE	Caps Assortment Kit			Ti, Nylon	1
335SBC	Smart Box Kit			Ti, Nylon	1
144AE	Analog for Lab			Stainless	2
044CAIN	Transfer for Impression			Stainless	2

COMPONENTS FOR TRI+ DIGITAL CAD/CAM

CAD/CAM Abutments - includes screw RS-TO



TO70-SCAN



TO70-07-F



TO70-07-NE



TO70-Blank-F

Product Code		Ø	GH	Material	Qty
TO70-SCAN	Scanbody for Laboratory and Intraoral Scanner			PEEK	1
TO70-07-F	Engaging Titanium Bonding-Base with TRI®-Friction	5.0 mm	0.7 mm	Ti-6Al-4V	1
TO70-07-NE	Non-Engaging Titanium Bonding-Base	5.0 mm	0.7 mm	Ti-6Al-4V	1
TO70-Blank-F	Titanium Milling Blank with TRI®-Friction			Ti-6Al-4V	1
RS-TO	Replacement Retaining Screw			Ti-6Al-4V	1

GH = Gingiva Height

Starting Drills

LD	Nr. 1	Lance Drill
RB-2.3	Nr. 2	Round Bur
DE	Nr. 3	Drill Extender

Surgical Drills

TPD-2.3S	Nr. 4	Pilot Drill - Short - 2.3 mm
TPD-2.3L	Nr. 5	Pilot Drill - Long - 2.3 mm
TND-3.3L	Nr. 6	TRI®-Narrow Final Drill - Long - 2.4/2.8 mm
TAP3.3	Nr. 7	TRI®-Narrow TAP
SBD-3.7S	Nr. 8	3.7 Soft Bone Drill - Short - 2.7/3.1 mmD
SBD-3.7L	Nr. 9	3.7 Soft Bone Drill - Long - 2.7/3.1 mmD
DBD-3.7L	Nr. 10	3.7 Dense Bone Drill - Long - 2.9/3.6 mmD
TAP3.7	Nr. 11	3.7 Tap
SBD-4.1S	Nr. 12	4.1 Soft Bone Drill - Short - 3.0/3.4 mmD
SBD-4.1L	Nr. 13	4.1 Soft Bone Drill - Long - 3.0/3.4 mmD
DBD-4.1L	Nr. 14	4.1 Dense Bone Drill - Long - 3.2/3.9 mmD
TAP4.1	Nr. 15	4.1 Tap
SBD-4.7S	Nr. 16	4.7 Soft Bone Drill - Short - 3.6/4.0 mmD
SBD-4.7L	Nr. 17	4.7 Soft Bone Drill - Long - 3.6/4.0 mmD
DBD-4.7L	Nr. 18	4.7 Dense Bone Drill - Long - 3.8/4.6 mmD
TAP4.7	Nr. 19	4.7 Tap
PP-L	Nr. 21	Paralell Pin - Long

Implant Drivers

VD-S	Classic Line TRI®-Vent Implant Driver - Short
VD-L	Classic Line TRI®-Vent Implant Driver - Long
ND-S	Classic Line TRI®-Narrow Implant Driver - Short
ND-L	Classic Line TRI®-Narrow Implant Driver - Long
OD-S	Classic Line TRI®-Octa Implant Driver - Short
OD-L	Classic Line TRI®-Octa Implant Driver - Long
MD-S-ISO	matrix® Implant Driver - Short
MD-L-ISO	matrix® Implant Driver - Long
MD-S-TW*	matrix® Implant Driver Short incl. Torque Wrench
MD-L-TW*	matrix® Implant Driver Long incl. Torque Wrench

Torque Wrench

TW	Torque Wrench
TW - Adapter	Torque Wrench Adapter

Prosthetic Tools

PD-L	Prosthetic Driver - Long
PD-S	Prosthetic Driver - Short
PD-ISO-L*	Prosthetic Driver for Motor - Long
PD-ISO-L-Ang*	Prosthetic Driver Angulated - Long
TORT	Removal Tool for Classic Line TRI®-Octa Friction Abutments
TRT	Removal Tool for Classic Line TRI®-Vent & Narrow Friction Abutments
487.ICE	TRI® Raptor Insertion and Extractor Tool for caps

Drill Stops


TDS6*	Drill Stop / 6 mm
TDS8*	Drill Stop / 8 mm
TDS10*	Drill Stop / 10 mm
TDS11*	Drill Stop / 11 mm
TDS13*	Drill Stop / 13 mm
TPDS*	Guided Surgery Sleeve for TPD-2.3L

Other


TSK	TRI® Surgical Kit - Complete
TSK-Tray	TRI® Surgical Kit - Tray

Prosthetic Kit

PK	TRI® Prosthetic Kit (includes: PD-ISO-L, PD-L, PD-S, TORT, TRT, TW, TW-Adapter)
PK-Tray	Prosthetic Kit-Tray
PD-ISO-L*	Prosthetic Driver for Motor - Long
PD-L	Prosthetic Driver - Long
PD-S	Prosthetic Driver - Short
TORT	Removal Tool for Classic Line TRI®-Octa Friction Abutments
TRT	Removal Tool for Classic Line TRI®-Vent & Narrow Friction Abutments
TW	Torque Wrench
TW - Adapter	Torque Wrench Adapter




Drill & Implant Lengths




Drill Sequence


3.3 mmD



3.75 mmD




4.1 mmD



4.7 mmD



Torque Requirements



TRI[®] Pod Handling

3 Options for Implant Insertion



Ratchet



Handpiece



By hand



Step 1
Open Lid



Step 2
Rotate To Remove
implant



Step 3
Remove Cover Screw

TO
TOUCH IT
IS TO
LOVE IT.

The implant packaging by TRI[®] integrates the latest engineering knowhow and innovation. Its revolutionary touchless delivery concept allows to pick-up the implant directly with the surgical hand-piece in one single step. Whilst guaranteeing the highest purity of the TRI[®] SBA surface. A true design innovation which will take your surgical performance to a new level. **Touch it, and you will love it!**

TRIACADEMY

To learn more about our products visit our TRI[®] Academy
Courses and Webinars:
academy@tri-implants.swiss




Become part of
the Community!



Join our
Facebook Group!

DESIGNED AND MANUFACTURED BY TRI[®] IN SWITZERLAND

Designed in the heart of Switzerland, all implants are produced in specialized manufacturing units for dental implants. Therefore benefiting from more than 25 years of medical technology production know-how - the only way to guarantee Swiss Precision with high quality.
Visit us in Switzerland: visit@tri-implants.swiss



TRI Dental Implants Int. AG
Switzerland
00800 3313 3313
www.tri.swiss



TRI® Dental Implants logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of TRI Dental Implants Int. AG.

TRI® Classic Line Catalog EN 2022-02/01



CE 0297 | FDA Registered