




TRI Dental Implants
Unique Products
May 2014

Through Research Innovative
www.tri-implants.com

TRI[®] - Vent & TRI[®] - Narrow

Unique products

- **The Healing Collar - Esthetic** has perforations in the upper range. It gives the surgeon the capability to thread in the suture material, to fix a vascular flap.



Height:

TRI-Vent & Narrow: 6mm


Diameter (Emergence Profil)

TRI-Vent
3,5mm
5mm

TRI-Narrow
3,5mm

TRI® - Vent & TRI® - Narrow

Unique products



Healing Collar – Esthetic

TVHC50-60-A



TVHC35-60-A

TNHC35-60-A

TRI[®] - Vent & TRI[®] - Narrow

Unique products

Healing Collar – Esthetic




Case by Dr. Spinato - Italy

TRI[®] - Vent & TRI[®] - Narrow

Unique products

Healing Collar – Esthetic



TRI[®] - Vent & TRI[®] - Narrow

Unique products

Surgical Cover Screw for Sinus Lift

Surgical Cover Screw
for Sinus Lift

5mm Durchmesser



Surgical Cover Screw -
standart

3.5mm Durchmesser



TRI[®] - Vent & TRI[®] - Narrow


Unique products

Surgical Cover Screw for Sinus Lift



TRI[®] - Vent & TRI[®] - Narrow

Surgical Cover Screw for Sinus Lift



TVSCS-SL
Surgical Cover Screw for Sinus Lift

TRI[®] - Vent & TRI[®] - Narrow

Unique products








Product Catalogue, page 2 & 3

SURGICAL COVER SCREW FOR SINUS LIFT

Catalog Number	Ø	Material	Qty	Sterile
TVSCS-SL	5.2 mm	Ti-6Al-4V	1	

Healing Collars, Implant Analogs

TISSUE HEALING – HEALING COLLAR

Catalog Number	Ø	Cuff Height	Material	Qty	Sterile
TVHC35-30 Healing Collar	3.5 mm	3.0 mm	Ti-6Al-4V	1	
TVHC35-45 Healing Collar	3.5 mm	4.5 mm	Ti-6Al-4V	1	
TVHC35-60 Healing Collar	3.5 mm	6.0 mm	Ti-6Al-4V	1	
TVHC35-60-A Healing Collar – Aesthetic	3.5 mm	6.0 mm	Ti-6Al-4V	1	
TVHC50-30 Healing Collar	5.0 mm	3.0 mm	Ti-6Al-4V	1	
TVHC50-45 Healing Collar	5.0 mm	4.5 mm	Ti-6Al-4V	1	
TVHC50-60 Healing Collar	5.0 mm	6.0 mm	Ti-6Al-4V	1	
TVHC50-60-A Healing Collar – Aesthetic	5.0 mm	6.0 mm	Ti-6Al-4V	1	